

Requisition for Sampling Containers and Supplies

CLIENT INFORMATION

Company		Contact Name							
Address									
Phone #		Email							
	Item cat #	Item/ Kit	# Bags	# Sponge	# Swab	# Bottles	Petri (2)	KIT 7211	
Food Sample Bags (usually 5 samples per pathogens) or 250g	6744	Sterile whirl Pak bag							
Sponge (usually 1 sponge per pathogen/site)	7888	Sterile sponge (10mL)							
Swab (usually 1 swab per site per pathogens)	4555	Sterile swab (10mL)							
Water microbiology (E. coli) bottles	2333	Sterile - Water bottle -250 ml							
Ice microbiology (E. coli) bottles	9654	Sterile – 2 bottles -250 ml							
Passive air sampling	3445	PCA and TSA							
Carcasses kit (include 3 sponges with Buffer peptone water and a plastic template)	7211	3 Sterile sponges with BPW							
Lab used only	Signature	Date	Sent Date						

Please choose the right bag or container or swab to be use for your analysis. Feel free to contact us if you have any questions. Call us at 705-586-1818 or contact us at info@capellainnovation.com